



Pressure Care Solutions



Prolevo is an innovative range of pressure redistribution and offloading devices for use in the prevention and management of pressure related injuries

About Prolevo

Prolevo is a division of TalarMade Limited, a company built by and for medical professionals with expertise in Orthotics, Podiatry, Diabetic Foot Care, and Physical Therapy. With over 40 years of knowledge and experience supplying quality products to the global healthcare market, we work closely with clinicians who have hands-on experience in practice and hospital settings.

By listening to and collaborating with experts in pressure redistribution, diabetic foot care, acute and community healthcare, and tissue viability, we develop products that provide **Pressure Care Solutions**, designed to support clinicians in their efforts to prevent pressure damage. We have a dedicated sales support team, alongside our clinical specialists, who can assist your organisation in developing an effective pressure ulcer prevention strategy.

Our mission is to create innovative products that are fit for purpose, with the goal of helping to reduce the risk of pressure ulcer development and support the healing process for patients with active pressure ulcers.

Our Range Includes:

• HeelPro	Page 3
• HeelPro Advance	Page 3
• FooSafe Prevention Boot	Page 4
• HeelSafe Pressure Redistribution Pad	Page 5
• SoleSafe Bed End Pressure Redistribution Pad	Page 5
• Safeguard Utility Pad	Page 6
• FootShield II Shoe	Page 7
• FootShield II Boot	Page 7
• FootShield II Wedge Shoe	Page 8
• FootShield Hex Insole	Page 8
• Ambulatory Pressure Relief Boot	Page 9
• Pressure Relief AFO	Page 9

Our helpful product selection pathways are on pages 10-11 if you need assistance selecting any products within our range.

To order our products or for more information, please contact us using the details below. Orders can be placed via telephone, email, NHSSC, or FP10, where appropriate based on the available ordering methods for each product.



01246 268456



info@talarmade.com



prolevo.com

HeelPro

For the prevention and management of heel pressure ulcers

- Offloads the heel
- Helps reduce incidence of heel pressure ulcers¹
- Heel aperture allows for easy inspection
- Universal sizing for easy application
- Vacuum packed to reduce storage space

* Supplied individually or in packs of 10

* Single patient use

Order Code	NHSSC Code	PIP Code	Qty
HPRO/U	GRD1502	409-1476	1
HPRO/U10	GAA123	---	10



Generous low friction cushioning for effective protection

HeelPro

Replacement Lycra Panels/Straps

- Can be used with both HeelPro and HeelPro Advance

Order Code	NHSSC Code	Qty
HPRO/RLP/2	---	2
HPRO/RLP	GRA1663	10



HeelPro

Stability Wedge

- Limits rotation of the leg
- Easily attaches to the boots
- Can be used with both HeelPro and HeelPro Advance

Order Code	NHSSC Code	Qty
HPRO/SW/U	GRA5182	1
HPRO/SW/U10	GRA5183	10

HeelPro Advance

- Offloads the heel
- Wipe clean cover which means it is ideal for use where there are existing wounds
- Heel aperture allows for easy inspection
- Universal sizing for easy application
- Vacuum packed to reduce storage space

* Supplied individually or in packs of 10

* Single patient use

Order Code	NHSSC Code	PIP Code	Qty
HPRO-AD/U	GAA092	414-6411	1
HPRO-AD/U10	GAA124	---	10



Ideal for use with existing wounds due to wipe clean cover

FootSafe Prevention Boot

- The FootSafe Boot may contribute to interventions in preventing and managing heel ulcers².
- Offloads the heel
- Supplied pre-inflated
- Anti-microbial material
- Easy close PU buttons and straps
- Durable

* Supplied in pairs or in packs of 10

* Multi patient use



Available in 3 sizes.

Style	Size	Qty (Pairs)	Product Code	NHSSC Code
Pre-Inflated	Small	1	PHPB-UC-PI-SM-P	GRA5228
Pre-Inflated	Small	5	PHPB-UC-PI-SM-P5	GRA5229
Pre-Inflated	Standard	1	PHPB-UC-PI-ST-P	GRA5230
Pre-Inflated	Standard	5	PHPB-UC-PI-ST-P5	GRA5231
Pre-Inflated	X-Large	1	PHPB-UC-PI-XL-P	GRA5232
Pre-Inflated	X-Large	5	PHPB-UC-PI-XL-P5	GRA5233

Sizing Guide

Size**	Shoe Size UK	Shoe Size EU
Small	2 - 7	34 - 40
Standard	7.5 - 11	41 - 45
X-Large	11.5+	46+

** If calf circumference is larger than normal (possibly due to oedema) consider fitting larger size than determined by shoe size.



Ideal for patients at risk of a fall

HeelSafe Pressure Redistribution Pad

- Redistributes Pressure
- Supplied inflatable or pre-inflated
- A solution for ambulatory patients and those at risk of falls when wearing boots
- An alternative choice for patient refusing or unsuitable for boots
- Multi-patient use
- Pressure mapped⁴
- Can be combined with SoleSafe
- One size - 80 x 66cm (W x L)

- * Supplied individually
- * Multi patient use

Style	Size	Pump	Qty	Size WxL (cm)	Product Code	NHSSC Codes
Covered	Pre-Inflated	---	Single	80 x 66	PHP-C-PI-OS	FES11326
Covered	Pre-Inflated	---	Pack 10	80 x 66	PHP-C-PI-OS10	FES11288
Covered	Inflatable	Yes	Single	80 x 66	PHP-C-IP-OS	FES11299
Covered	Inflatable	Yes	Pack 10	80 x 66	PHP-C-IP-OS10	FES11253
Covered	Inflatable	No	Single	80 x 66	PHP-C-NP-OS	FES11323
Covered	Inflatable	No	Pack 10	80 x 66	PHP-C-NP-OS10	FES11242

SoleSafe Bed End Pressure Redistribution Pad


- Redistributes Pressure
- Supplied inflatable or pre-inflated
- Protects the soles of feet
- Allows objects hung on the end of the bed to be concealed
- Pressure mapped³
- Can be combined with HeelSafe above
- One size - 80 x 28cm (W x L)

- * Supplied individually
- * Multi patient use



Ideal for patients who may slide down the bed

Style	Size	Pump	Qty	Size WxL (cm)	Product Code	NHSSC Codes
Covered	Pre-Inflated	---	Single	80 x 28	PSP-C-PI-OS	FES11244
Covered	Pre-Inflated	---	Pack 10	80 x 28	PSP-C-PI-OS10	FES11266
Covered	Inflatable	Yes	Single	80 x 28	PSP-C-IP-OS	FES11243
Covered	Inflatable	Yes	Pack 10	80 x 28	PSP-C-IP-OS10	FES11289
Covered	Inflatable	No	Single	80 x 28	PSP-C-NP-OS	FES11291
Covered	Inflatable	No	Pack 10	80 x 28	PSP-C-NP-OS10	FES11295



L 50 x W 30 x H 2.5 cm
Product can be cut to shape



Safeguard Utility Pad

The Safeguard Utility Pad is designed to assist in the prevention of pressure ulcers as part of a patient's overall pressure ulcer prevention care program.

The Utility Pad redistributes pressure and can be used on bony prominences such as knees, elbows, spinous processes and heels.

- * Supplied in pairs
- * Single patient use



Order Code	PIP Code	NHSSC Code	Qty
SGFUP/U	413-2544	FES11251	1
SGFU/U10	---	FES11308	10

Footshield II Shoe

- Evidence based design⁵
- Left and right models to maximise function and outcomes
- Reduces forefoot loads and pressure
- Accommodates bulky dressings and swollen feet
- Comes as standard with memory foam insole (optional Hex offloading insole on page 8)
- Removable toe cap



* Supplied individually

Size	UK Shoe Size	Length (mm)	LEFT			RIGHT		
			Product Code	NHSSC Code	PIP Code	Product Code	NHSSC Code	PIP Code
X-Small	1.5 - 3	240	FSS-II/XSL	ERP5998	---	FSS-II/XSR	ERP5999	---
Small	3.5 - 5	257	FSS-II/SL	ERP5994	422-5991	FSS-II/SR	ERP5995	422-6296
Medium	5.5 - 7	273	FSS-II/ML	ERP5992	422-6478	FSS-II/MR	ERP5993	422-6502
Large	7.5 - 9	290	FSS-II/LL	ERP5990	422-6270	FSS-II/LR	ERP5991	422-6288
X-Large	9 - 11	307	FSS-II/XLL	ERP5996	422-6486	FSS-II/XLR	ERP5997	422-6494
WITH HEX INSOLE								
X-Small	1.5 - 3	240	FSS-IIH/XSL	ERP5988	---	FSS-II/XSR	ERP5989	---
Small	3.5 - 5	257	FSS-IIH/SL	ERP5984	---	FSS-IIH/SR	ERP5985	---
Medium	5.5 - 7	273	FSS-IIH/ML	ERP5982	---	FSS-IIH/MR	ERP5983	---
Large	7.5 - 9	290	FSS-IIH/LL	ERP5980	---	FSS-IIH/LR	ERP5981	---
X-Large	9 - 11	307	FSS-IIHH/XLL	ERP5986	---	FSS-IIH/XLR	ERP5987	---

Removable toe cap

Scan QR to view our Clinical Evidence Brochure



Footshield II Boot

- Features as above for the Footshield II Shoe - but with an extended heel and greater height

* Supplied individually

Size	UK Shoe Size	Length (mm)	LEFT		RIGHT	
			Product Code	NHSSC Code	Product Code	NHSSC Code
X-Small	1.5 - 3	240	FSB-II/XSL	ERP6018	FSB-II/XSR	ERP6019
Small	3.5 - 5	257	FSB-II/SL	ERP6014	FSB-II/SR	ERP6015
Medium	5.5 - 7	273	FSB-II/ML	ERP6012	FSB-II/MR	ERP6013
Large	7.5 - 9	290	FSB-II/LL	ERP6010	FSB-II/LR	ERP6011
X-Large	9 - 11	307	FSB-II/XLL	ERP6016	FSB-II/XLR	ERP6017
WITH HEX INSOLE						
X-Small	1.5 - 3	240	FSB-IIH/XSL	ERP5978	FSB-II/XSR	ERP5979
Small	3.5 - 5	257	FSB-IIH/SL	ERP6034	FSB-IIH/SR	ERP5975
Medium	5.5 - 7	273	FSB-IIH/ML	ERP6032	FSB-IIH/MR	ERP6033
Large	7.5 - 9	290	FSB-IIH/LL	ERP6030	FSB-IIH/LR	ERP6031
X-Large	9 - 11	307	FSB-IIH/XLL	ERP5976	FSB-IIH/XLR	ERP5977



Accommodates bulky dressings



Removeable memory foam insole

Footshield II Wedge Shoe

- Forefoot off-loader
- Removable memory foam insole
- Accommodates bulky dressings
- Can be used with optional Hex offloading insole below
- Removable toe cap


* Supplied individually

Size	UK Shoe Size	Length (mm)	LEFT		RIGHT	
			Product Code	NHSSC Code	Product Code	NHSSC Code
X-Small	1.5 - 3	240	FSWS-II/XSL	ERP6028	FSWS-II/XSR	ERP6029
Small	3.5 - 5	257	FSWS-II/SL	ERP6024	FSWS-II/SR	ERP6025
Medium	5.5 - 7	273	FSWS-II/ML	ERP6022	FSWS-II/MR	ERP6023
Large	7.5 - 9	290	FSWS-II/LL	ERP6020	FSWS-II/LR	ERP6021
X-Large	9 - 11	307	FSWS-II/XLL	ERP6026	FSWS-II/XLR	ERP6027
WITH HEX INSOLE						
X-Small	1.5 - 3	240	FSWS-IIH/XSL	ERP5988	FSWS-IIH/XSR	ERP5989
Small	3.5 - 5	257	FSWS-IIH/SL	ERP5984	FSWS-IIH/SR	ERP5985
Medium	5.5 - 7	273	FSWS-IIH/ML	ERP5982	FSWS-IIH/MR	ERP5983
Large	7.5 - 9	290	FSWS-IIH/LL	ERP5980	FSWS-IIH/LR	ERP5981
X-Large	9 - 11	307	FSWS-IIH/XLL	ERP5986	FSWS-IIH/XLR	ERP5987

Footshield Hex Insole

- Removable hexagonal pegs create a void to offload ulcerated or at risk sites
- Cushioned polyurethane top cover

* Supplied individually



Size	UK Shoe Size	Length (mm)	SHOE/BOOT		WEDGE SHOE	
			Product Code	NHSSC Code	Product Code	NHSSC Code
X-Small	1.5 - 3	240	FSHISB/XS	ERF616	FSHIWS/XS	ERF614
Small	3.5 - 5	257	FSHISB/S	---	FSHIWS/S	ERF612
Medium	5.5 - 7	273	FSHISB/M	---	FSHIWS/M	ERF611
Large	7.5 - 9	290	FSHISB/L	---	FSHIWS/L	ERF610
X-Large	9 - 11	307	FSHISB/XL	ERF615	FSHIWS/XL	ERF613

Ambulatory Pressure Relief Boot

- Offloads pressure from heel
- Soft, breathable and washable liner
- Adjustable length toe support
- Anti-rotation bar
- Available in two sizes

* Supplied individually

Size	Calf Circ. (cm)	Order Code	NHSSC Code
Standard	Max. 45	APRB/ST	GRD3625
Large	38 - 55	APRB/L	GRD3624

Sole plate can easily be removed and reapplied



Pressure Relief AFO

- Universal size
- Offloads pressure from heel
- Tri-laminate material allows air to circulate and helps wick moisture
- Adjustable length toe support
- Padded ankle strap reinforces foot position
- Supplied with spare liner set
- Non-slip sole allows for ambulation

* Supplied individually

Size	Product Code	NHSSC Code
Universal	AUPRNS/U	GRD3627

Hip rotation control bar



AT RISK PATIENTS
Product choice pathway

AMBULANT PATIENT

NON-AMBULANT PATIENT

Fit device to bed

Patient worn or fit to bed



ACTIVE FOOT OR HEEL ULCER PATIENTS
Product choice pathway

AMBULANT PATIENT

NON-AMBULANT PATIENT

Rigid frame device

Non-rigid device



How to Order:

To order our products or for more information, please contact us using the details below.

Orders can be placed via telephone, email, NHSSC, or FP10, where appropriate based on the available ordering methods for each product.

 01246 268456

 info@talarmade.com

 prolevo.com



prolevo

Pressure Care Solutions



TalarMade
Made with clinical excellence in mind

Prolevo is a division of TalarMade Ltd
Edinburgh House
Millennium Way
Dunston Technology Park
Chesterfield
S41 8ND



01246 268456



info@talarmade.com



prolevo.com

References

1. Reducing PU incident rates by giving the run around. Hayley Jones, Tissue Viability Advanced Nurse Specialist & Gemma McGrath, Tissue Viability Clinical Nurse Specialist, Medway NHS Foundation Trust. Poster, Wounds UK Annual Conference, November 2017.
2. Evaluation of a pressure ulcer prevention offloading boot in the acute setting. Deborah Gleeson, Lead Nurse, Tissue Viability. August 2017.
3. Evaluation of Prolevo SoleSafe. Duncan Bain Consulting, available on request.
4. Evaluation of Prolevo HeelSafe. Duncan Bain Consulting, available on request.
5. Effect of rocker sole design features on forefoot plantar pressures in people with and without diabetes. Chapman et al, Clinical Biomechanics, 28(6) 679-85, July 2013.

ART-0124-02